

PLEASE NOTE

The dimensions highlighted in "gold" are the common ordering measurement



**PESSARY REMOVER**

These devices are useful in removing pessaries, and are only available to physicians and nurse practitioners.

**BIOTEQUE'S PESSARY FITTING SET**

An aid in selecting the proper patient pessary.



**Ring**



SIZE	INCH	MM
1	2.00	51
2	2.25	57
3	2.50	64
4	2.75	70
5	3.00	76
6	3.25	83
7	3.50	89
8	3.75	95
9	4.00	102

w & w/o support

**Donut**



SIZE	INCH	MM
0	2.00	51
1	2.25	57
2	2.50	64
3	2.75	70
4	3.00	76
5	3.25	83
6	3.50	89

**Gelhorn**



SIZE	INCH	MM
0	1.50	38
1	1.75	44
2	2.00	51
3	2.25	57
4	2.50	64
5	2.75	70
6	3.00	76
7	3.25	83
8	3.50	89

w drains

**Short Stem Gelhorn**



SIZE	INCH	MM
0	1.50	38
1	1.75	44
2	2.00	51
3	2.25	57
4	2.50	64
5	2.75	70
6	3.00	76
7	3.25	83
8	3.50	89
9	3.75	95

w drains

**Shatz**



SIZE	INCH	MM
0	1.50	38
1	1.75	44
2	2.00	51
3	2.25	57
4	2.50	64
5	2.75	70
6	3.00	76
7	3.25	83
8	3.50	89

**Marland**



SIZE	INCH	MM
2	2.25	57
3	2.50	64
4	2.75	70
5	3.00	76
6	3.25	83
7	3.50	89
8	3.75	95

w & w/o support

**Ring with Knob**



SIZE	INCH	MM
1	2.00	51
2	2.25	57
3	2.50	64
4	2.75	70
5	3.00	76
6	3.25	83
7	3.50	89

w & w/o support

**Dish**



SIZE	INCH	MM
0	1.97	50
1	2.17	55
2	2.36	60
3	2.56	65
4	2.75	70
5	2.95	75
6	3.15	80
7	3.35	85

w & w/o support

**Oval**



SIZE	INCH	MM
1	2.00	51
2	2.25	57
3	2.50	64
4	2.75	70
5	3.00	76
6	3.25	83
7	3.50	89
8	3.75	95
9	4.00	102

w & w/o support

**Cube**



SIZE	INCH	MM
0	0.98	25
1	1.14	29
2	1.30	33
3	1.46	37
4	1.61	41
5	1.77	45
6	1.97	50
7	2.20	56
8	2.48	63
9	2.75	70
10	2.95	75

w & w/o drains

**Gehring**



SIZE	INCH	MM
0	1.97	50
1	2.17	55
2	2.36	60
3	2.50	65
4	2.75	70
5	2.95	75
6	3.15	80
7	3.35	85
8	3.54	90

**Cup**



SIZE	INCH	MM
0	1.97	50
1	2.17	55
2	2.36	60
3	2.56	65
4	2.75	70
5	2.95	75
6	3.15	80
7	3.35	85

w & w/o support

**Hodge**



SIZE	INCH	MM
0	2.60	65
1	2.75	70
2	2.95	75
3	3.15	80
4	3.37	85
5	3.54	90
6	3.74	95
7	3.94	100
8	4.14	105
9	4.33	110

w & w/o support

**Inflatable Donut**



SIZE	INCH	MM
2	2.00	51
3	2.25	57
4	2.50	64
5	2.75	70

**Flexi-Shelf**



SIZE	INCH	MM
6	2.25	57
8	2.50	64
10	2.75	70
12	3.00	76

**Incontinence Ring**



SIZE	INCH	MM
0	1.75	44
1	2.00	51
2	2.25	57
3	2.50	64